

T-Link mobile phone neck loop



The T-Link loop is placed around the neck. The other end is connected to the hands free socket of the mobile phone. It has a built-in microphone for hands free operation. When using the phone you switch their hearing aid/s to the 'T' position to receive the incoming speech. T-Links are available with a wide range of connectors and are therefore suitable for most mobile phones.

Supplier [Connevans](#)

Manufacturer

Product Code C01

Sensory Features Hands free system that can be used with hearing aids

Price £63.00 (incl. VAT)

Order Online This product can be ordered from [Deaf Equipment](#)

Catalogue A catalogue can be ordered from:

Connevans Ltd
54 Albert Road North
Reigate, Surrey
RH2 9YR

<http://www.deafequipment.co.uk>

More Information For more information, please contact askaquestion@connevans.com

©2012 Sense

Sense, 101 Pentonville Road, London N1 9LG, Tel.: 0845 127 0060, Text 0845 127 0062, email: info@sense.org.uk

Charity No. 289868

Built on [Cubik](#)

[Sign In](#) [Sign Up](#)